

Artic pipeline

- Description
- Setup
 - Set up Guppy
 - Set up conda environment
- Usage
 - Input sample sheet
 - Input reads
- Devel

Description

...

Setup

Set up Guppy

Download the appropriate version of guppy from [Oxford Nanopore](#) (requires registration, which is free), e.g. `ont-guppy_6.4.2_linux64.tar.gz` (GPU) or `ont-guppy-cpu_6.4.2_linux64.tar.gz` (CPU).

Or get it from

```
# CPU
```

```
wget https://mirror.oxfordnanoportal.com/software/analysis/ont-guppy-cpu_6.4.2_linux64.tar.gz
```

```
# GPU
```

```
wget https://mirror.oxfordnanoportal.com/software/analysis/ont-guppy-gpu_6.4.2_linux64.tar.gz
```

Extract files:

```
tar zxvf ont-guppy_6.4.2_linux64.tar.gz
```

Then add the `bin` directory to your `PATH` variable:

```
export PATH=/full/path/to/ont-guppy_6.4.2_linux64/bin:$PATH
```

To permanently have guppy available on your `PATH`, add the command above to the file `~/.bashrc`.

If you don't or you can't edit your `PATH`, use option `--guppy-path` in `artic-smk.py` to point to the guppy `bin` directory. E.g. `--guppy-path /path/to/ont-guppy_6.4.2_linux64/bin`

Set up conda environment

- Install [conda](#), [mamba](#), and configure for [bioconda](#).
- Create a dedicated environment for this pipeline

```
conda create --yes -n artic-smk
```

```
conda activate artic-smk
```

```
mamba install --yes --file requirements.txt -n artic-smk
```

Usage

The following command should work as-is using the test data. It will process the given `fast5` directory according to `sample_sheet.tsv`. Since option `--dry-run` is set it will only print what would be executed, remove it for the real processing.

```
./artic-smk.py --sample-sheet test/data/sample_sheet.tsv \
  --fast5-dir test/data/fast5 \
  --genome-name my-genome \
  --output test_out \
  --dry-run
```

Run `./artic-smk.py -h` to see the list of available options. The following printout may be out of date:

Run artic pipeline

optional arguments:

```
-h, --help                show this help message and exit
--sample-sheet SAMPLE_SHEET, -s SAMPLE_SHEET
                           Tab-separated file of samples, barcodes, and
                           other sample-specific options. See online
                           documentation for details [required]
--output OUTPUT, -o OUTPUT Output directory [artic-out]
--medaka-scheme-directory MEDAKA_SCHEME_DIRECTORY, -sd MEDAKA_SCHEME_DIRECTORY
                           Path to scheme directory [primer-schemes]
--fast5-dir FAST5_DIR, -f5 FAST5_DIR
                           Directory of fast5 file. Typically the
                           Nanopore run directory
--fastq-dir FASTQ_DIR, -fq FASTQ_DIR
                           Input alternative to fastq5-dir: Directory
                           of demultiplexed fastq files. fastq-dir
                           contains subdirectories named after the
                           sample barcodes and containing the
                           respective fastq files
--genome-name GENOME_NAME, -g GENOME_NAME
                           Name for consensus genome [genome]
--guppy-config GUPPY_CONFIG
                           For fast5 input: Configuration for
                           guppy_basecaller
                           [dna_r9.4.1_450bps_fast.cfg]
--guppy-barcode-kit GUPPY_BARCODE_KIT
                           For fast5 input: Barcode kit passed to
                           guppy_barcode [EXP-NBD104]
--guppy-basecaller-opts GUPPY_BASECALLER_OPTS
                           Additional options passed to
                           guppy_basecaller as a string with leading
                           space e.g. " --num_caller 10" []
--guppy-path GUPPY_PATH    Full path to guppy bin directory. Leave
                           empty if guppy is already on your search
                           PATH []
--min-length MIN_LENGTH, -L MIN_LENGTH
                           Ignore reads less than min-length [350]
--medaka-model MEDAKA_MODEL
                           Model to use for medaka [r941_min_fast_g303]
--jobs JOBS, -j JOBS      Number of jobs to run in parallel [1]
--dry-run, -n             Run pipeline dry-run mode
--snakemake-opts SNAKEMAKE_OPTS, -smk SNAKEMAKE_OPTS
                           Additional options to snakemake as a string
                           with leading space e.g. " --rerun-incomplete
                           -k" []
--version, -v             show program's version number and exit
```

Input sample sheet

This is a tabular file tab or comma separated with first non-skipped line as header. Lines starting with '#' are skipped. Columns are:

Column	Description
sample	Sample name. Avoid names with spaces or special characters (dots, underscores, hyphens are ok)
barcode	Sample barcode

Additional columns are ignored

Input reads

- **Options 1** A directory of **fast5** files that will be passed to **guppy_basecaller** and **guppy_barcode**. Typically this is the output of the Nanopore run. Use **--fast5-dir/-f5** option to start from here.
- **Options 2** A directory of **fastq** files already demultiplexed and ready for further processing. Use **--fastq/fq** option to start from here, guppy installation is not required. Fastq-dir contains subdirectories named after the sample barcodes. They don't need to be real barcode names as long as they match the sample sheet column **barcode**. Each subdirectory can contain multiple fastq files, possibly gzip'd. This is the test data example:

```
test/data/fastq/  
  barcode01  
    tvla1_run2.fastq.gz  
  barcode02  
    tvla1_run1.fastq.gz  
  barcode04  
    dummy04.fastq.gz  
  barcode06  
    dummy06.fastq.gz
```

Devel

To run test suite:

```
./test/test.py
```

Compile this markdown to pdf:

```
pandoc -V colorlinks=true -V geometry:margin=1in README.md -o README.pdf
```